


PerfectSeal[®]
Healthcare



PerfectSeal[®] | DENSUN

Add: A2-1 OF # 5 PLANT , NO.96 MUJIN ROAD, NINGBO
HIGH-TECH ZONE, NINGBO CITY, ZHEJIANG PROVINCE.

Tel: +86-574-8756 4800

Fax: +86-574-8756 4717

E-mail: sales@perfectseal.com

Http: www.perfectseal.com



2023-2024
MEDICAL EQUIPMENT



国家高新技术企业

COMPANY PROFILE



Zhejiang PerfectSeal New Material Technology Co., Ltd. specializes in the design, development, manufacturing, and sales of professional Medical device products. The company's main products include syringes, infusion sets, injection needles, insulin needles, disposable wound dressings and first-aid products. The company has a modern 100,000-grade purification workshop in line with GMP standards, as well as perfect R & D and quality control system, advanced inspection machine, ISO13485 quality system management certification, to ensure that our products meet the relevant regulatory requirements, to meet the strict requirements of customers in different regions and countries.

The company has long insisted on the principle of quality first and quality first and strengthened scientific management. The company continues to develop and grow, and sincerely welcomes all major medical institutions and companies in China to sincerely cooperate with us and make their own outstanding contributions to the cause of human health.



Quality Management:

ISO 13485 standard is an effective solution to meet the comprehensive requirements of QMS. The adoption of ISO 13485 is an effective foundation for manufacturers to meet the requirements of medical device regulations and their commitment to the safety and quality of their products. It ensures that medical devices are designed, developed, manufactured, installed and delivered in accordance with the requirements on an ongoing basis and that they are safe for use as intended.



Corporate Philosophy:

The company insists on the purpose of "quality first, sincere service, care for health", adhering to the vision of science and technology to make life better, and long-term focus on the health care industry. PerfectSeal is willing to work with domestic and foreign counterparts to make its own outstanding contributions to the further development of medical and health care in China and the world.

Co-existence, win-win and common development" is the constant tenet of PERFECTSEAL

PerfectSeal's Development History



Contents

Disposable Syringe

Disposable Syringe.....	05
Insulin Syringe	05
Destruct Syringe	05

Disposable Syringe Needle

I.V. Cannula.....	11
Sterile Irrigation Needle	11
Beauty Needle	11
Insulin Pen Needle	11

Disposable Infusion set

Disposable Infusion Set	17
Disposable Infusion Set ("Y" Injection port)	17
Disposable Infusion Set (Precise).....	17

Respiratory Products

Disposable Anaesthesia Mask	21
Disposable Nebulizer Mask	21
Disposable Oxygen Mask	21

Endotracheal Tubes

Endotracheal Tubes.....	24
Endotracheal Tubes.....	24

Urinary Drainage Bag

Urinary Drainage Bag.....	24
---------------------------	----

Wound Care Products

Hydrocolloid Wound Dressing.....	25
Foam Dressing.....	25
Film plus.....	26
PORE/Self-Adhesive Wound Dressing.....	27
Strips-Wound skin closures.....	27
Sports tape - kinesiology tape.....	28
Silicone Foam Dressing.....	28

First Aid Bag

First Aid Bag.....	29
--------------------	----

Pet First Aid Bag

Pet First Aid Bag.....	33
------------------------	----

Production Scene



1 Production Scene



2 Production Scene




3 Production Scene



4 Production Scene



5 Production Scene



6 Production Scene



1 (Luer Slip)



2 (Luer Lock)



4



5



6

Disposable Syringes

Disposable Syringe

Part No.:

Middle head type and offset head type, straight mouth type and screw type, two and three piece type; Soft medium package, hard medium package; With or without a needle.

◆ Product specifications

From 1ml to 60ml


◆ Syringe with needle Injection needle specification:

Needle size from 0.45mm to 1.2mm

Part No.	Capacity	Needle	Application	Part No. (International)	Part No. (China)
PDS-001	1ml	0.45mm	Subcutaneous injection	26G	No.4.5
PDS-002	2ml	0.5/0.6mm	Venous blood collection	26G/24G	No.5/No.6
PDS-003	3ml	0.5/0.6mm	Venous blood collection	26G/24G	No.5/No.6
PDS-004	5ml	0.6/0.7mm	Venous blood collection	24G/22G	No.6/No.7
PDS-005	10ml	0.8/1.2mm	Aspirate the liquid	18G	No.8/No.12
PDS-006	20ml	1.2/1.6mm	Aspirate the liquid	18G	No.12/No.16
PDS-007	30ml	1.2mm	Aspirate the liquid	18G	No.12
PDS-008	50ml	1.2mm	Aspirate the liquid	18G	No.12/No.16
PDS-009	60ml	1.2mm	Aspirate the liquid	18G	No.12/No.16

Constitute

The structure drawing is as follows:



Attention

- ◆ This product and syringe are used for injection or extraction of liquid medicine in human skin, skin, muscle, vein, etc.
- ◆ This product is sterilized by ethylene oxide, sterile, no heat source;
- ◆ Open the package, please use immediately to avoid contamination;
- ◆ Single packaging damage, protective cover off, abnormal needle tip, strictly prohibited use;
- ◆ Use only once, destroy after use;
- ◆ The needle tube should not be bent for use, otherwise there is a risk of needle breakage;

Application



Disposable Syringe



Needle	Capacity	Brown
26G(0.45)×13mm	1ml	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/> Luer Lock <input type="checkbox"/> Luer Slip		

Needle	Capacity	Brown
26G(0.45)×13mm	1ml	<input type="checkbox"/>
<input type="checkbox"/> Luer Lock <input checked="" type="checkbox"/> Luer Slip		

Disposable Syringe



Needle	Capacity	Orange
25G(0.50)×25mm	3ml	
<input checked="" type="checkbox"/> Luer Lock <input type="checkbox"/> Luer Slip		



Needle	Capacity	Black
22G(0.60)×38mm	5ml	
<input type="checkbox"/> Luer Lock <input checked="" type="checkbox"/> Luer Slip		



Needle	Capacity	Light green
21G(0.80)×38mm	10ml	
<input type="checkbox"/> Luer Lock <input checked="" type="checkbox"/> Luer Slip		



Needle	Capacity	Pink
18G(1.20)×38mm	20ml	
<input checked="" type="checkbox"/> Luer Lock <input type="checkbox"/> Luer Slip		

Insulin Syringe

Product features:

- ◆ Ultra micro needle, reduce injection pain.



Part No.	Capacity
PIS-001	0.3ml, 0.5ml, 1.0ml



Part No.	Capacity
PIS-001	0.3ml, 0.5ml, 1.0ml

Destruct Syringe

Product features:

- ◆ Self-destruct function automatically starts after injection, effectively prevent secondary use. There is a clamp design at the entrance of the jacket. After injection, the clamp will automatically resist the core rod, effectively preventing the core rod from being pulled out of the jacket and used again (the standard requires 100N pulling force), realizing the self-destruction function and avoiding cross infection.
- ◆ Simple product structure design, effective control of drug residue.



Part No.	Capacity
PSS-001	0.1ml, 0.2ml, 0.5ml, 1ml, 2ml, 3ml, 4ml, 5ml



Part No.	Capacity
PSS-002	0.5ml, 1ml, 2ml, 3ml, 5ml, 10ml

Syringe Needle




Syringe Needle

Product model:
From 0.45 mm to 1.2 mm

Product features


- ◆ High quality austenitic SUS304 stainless steel; • The syringe adopts the world's leading thin-wall design, with large inner diameter and large flow rate;
- ◆ Needle tip design with excellent sharpness, rapid insertion, reduce pain, minimal tissue damage;
- ◆ 6:100 conical connector manufactured in accordance with international standards to ensure accurate size and good compatibility with other medical equipment;
- ◆ Needle seat color recognition standard, easy to distinguish the use. Translucent design, blood reflux at a glance;
- ◆ The protective shell can easily and effectively shield the used needle after use, so as to ensure that the needle will not accidentally stab the medical staff, patients and other people, so as to achieve safety and reduce infection.



Top crops



Elbow



Needle without cut



Round needle

International Reference Table

Specification		Color code	Length metric system (mm)							Wall thickness		Edge angle	
Metric system	Gauge		13	16	19	25	32	38	50	RW	TW	LB	SB
0.25	31G	Blue											
0.30	30G	Yellow											
0.33	29G	Red											
0.36	28G	Blue-green											
0.40	27G	Medium grey											
0.45	26G	Brown											
0.50	25G	Orange											
0.55	24G	Medium purple											
0.60	23G	Dark blue											
0.70	22G	Black											
0.80	21G	Dack green											
0.90	20G	Yellow											
1.10	19G	Creamy											
1.20	18G	Pink											

Notes: In addition to the above common specifications, special specifications can be customized according to customer requirements.

Disposable Syringe Needle



Item No.	PSN-001	PSN-002	PSN-003	PSN-004	PSN-005
SIZE	14g-26g	14g-26g	14g-26g	14g-26g	14g-26g

Disposable Syringe Needle

Product Features:

- Special tapered & kink resistant FEP radiopaque catheter for easy insertion with optimal flow rate
- Special flexible wings to provide easy & secure fixation
- Manufactured from tested Biocompatible materials offering longer indwelling time
- Rounded grip on the needle hub for better support during insertion.
- Injection port with non-return silicon valve for intermittent medication
- Transparent flashback chamber for easy visual medication



Item No.	Size	Cs (PCS)
PIV-001	14g,16g,18g,20g,22g,24g,26g	1 pcs/bag
PIV-002	14g,16g,18g,20g,22g,24g,26g	1 pcs/bag

Disposable Sterile Irrigation Needle

Product Features

- Made from high-quality austenitic stainless steel;
- The needle tube adopts the world-leading thin-wall design with large inner diameter and large flow rate;
- Needle base color identification specification for easier differentiation;
- The 6:100 threaded and non-threaded conical connector manufactured according to the international standard ensures accurate size and good compatibility with other medical devices;
- The needle tube's bending size, bending angle, needle tip shape, etc. can be customized according to customer requirements;



Item No.	Size	Cs (PCS)
PIN-001	14g,16g,18g,20g,22g,24g,26g	1 pcs/bag

Disposable Beauty Needle

Product Features

- High quality austenitic stainless steel;
- The syringe adopts the world's leading thin-wall design, with large inner diameter and large flow;
- Needle seat color recognition specification, easy to distinguish;
- 6:100 threaded and non-threaded tapered connectors manufactured according to international standards to ensure accurate size and good compatibility with other medical equipment;
- The bending size of needle tube, bending angle, shape of needle tip can be customized according to customer requirements.



PBN-001



Item No.	Size	Cs (PCS)
PBN-001	14g,16g,18g,20g,22g,24g,26g	1pcs/bag

Disposable Insulin Pen Needle

Product Features:

- This product consists of sealing strip, needle seat, needle tube, small sheath and large sheath, for use in supporting the insulin pen on subcutaneous injection in the human body;
- The product meets the requirements of ISO 1168. 2 standard, having passed biological evaluation test and the compatibility test with insulin solution;
- The needle tip of the product is particularly tiny and sharp, causing mild pain in injection, little skin damage, and minimal residual liquid;
- With complete specifications, it is suitable for supporting various insulin pens on-the-shelf in the international market,



PPN-001

Specifications		Color code	Needle tube length Metric system(mm)			
Metric system	Gauge		12	8	5	5
0.20	33G	Purple				
0.23	32G	Green				
0.25	31G	Blue				
0.30	30G	Yellow				
0.33	29G	Pink				

Item No.	Size	Cs (PCS)
PPN-001	14g,16g,18g,20g,22g,24g,26g	1pcs/bag

Central Venous Catheter

Product Features:

- Movable clamp allows anchorage at puncture site regardless of depth of catheter, which minimizes trauma and irritation to the puncturing site. Depth marking assists in accurate placement of central venous catheter from the right or left subclavian or jugular vein. Soft tip reduces trauma to vessel, minimizing vessel erosion, hemothorax and cardiac tamponade. Single, double, triple and quad lumen is available for choice.

Item No.	Size	Cs (PCS)
Central Venous Catheter	Standard kits/Optional Compound Kits	1pcs/bag

Standard kits include:

- Central Venous Catheter
- Guide-wire
- Vessel Dilator
- Clamp
- Fastener: Catheter Clamp
- Introducer Needle
- Introducer Syringe
- Injection Needle
- Injection Cap

Optional Compound Kits include:


- Central Venous Catheter Standard Kit Accessories
- 5ml Syringe
- Surgical Gloves
- Clamp
- Surgical Pledget
- Surgery Sheet
- Surgery Towel . Sterile Brush
- Gauze Pad
- Suture of Needle
- Wound Dressing
- Scalpel



Disposable Infusion Set

Item No.	Size	Cs (PCS)
PII-001	150cms, Luer lock with Gauge 21G Needle	1pcs/bag
PII-002	150cms, Luer Slip with Gauge 21G Needle	1pcs/bag
PII-003	150cms, Luer lock with Gauge 21G Needle	1pcs/bag
PII-004	150cms, Luer lock with Gauge 21G Needle	1pcs/bag

Constitute




Disposable Infusion Set (Luer lock)

Product Features:

- Transparent & flexible standard drip chamber
- 15 micro fluid filter in drip chamber
- Sharp and easy piercing air vented spike
- Bacteria retentive air inlet with snap on cap
- Approximately 20 drops/ml
- Also available without airvent

Soft and kink resistant PVC tubing ensures uniform flow rate
 Smooth roller clamp facilitates easy, & controlled adjustment of fluid rates
 Standard Tube: 150cms
 Luer lock with Gauge 21G Needle



Disposable Infusion Set

Disposable Infusion Set (Luer Slip)

Product Features:

- Transparent & flexible standard drip chamber
- 15 micro fluid filter in drip chamber
- Sharp and easy piercing air vented spike
- Bacteria retentive air inlet with snap on cap
- Approximately 20 drops/ml
- Also available without airvent

Soft and kink resistant PVC tubing ensures uniform flow rate
 Smooth roller clamp facilitates easy, & controlled adjustment of fluid rates
 Standard Tube:150cms
 Luer Slip with Gauge 21G Needle



PII-002

Disposable Infusion Set ("Y" Injection port)

Product Features:

- Transparent & flexible standard drip chamber
- 15 micro fluid filter in drip chamber
- Sharp and easy piercing air vented spike
- Bacteria retentive air inlet with snap on cap
- Approximately 20 drops/ml
- Also available without airvent

Soft and kink resistant PVC tubing ensures uniform flow rate
 Smooth roller clamp facilitates easy, & controlled adjustment of fluid rates
 Standard Tube :150cms
 Luer lock "Y" Injection port (Latex Free) designed for intermittent medications



PII-003

Disposable Infusion Set (Precise)

Product Features:

- Transparent & flexible standard drip chamber
- 15 micro fluid filter in drip chamber
- Sharp and easy piercing air vented spike
- Bacteria retentive air inlet with snap on cap
- Approximately 20 drops/ml
- Also available without airvent

Soft and kink resistant PVC tubing ensures uniform flow rate
 Smooth roller clamp facilitates easy, & controlled adjustment of fluid rates
 Standard Tube :150cms
 Luer lock "Y" Injection port (Latex Free) designed for intermittent medications



PII-004

Hemodialysis Blood Set

Product Features:

- All tubes are made from medical grade, and all components manufactured in original.
- Pump tube: With high elasticity and medical grade PVC, the shape of tube remains the same after continuous pressing of 10 hours.
- Drip Chamber: several sizes of drip chamber available.
- Dialysis Connector: Extra large designed dialyzer connector is easy to operate.
- Clamp: Clamp is made of hard poastic and designed large and thicker to guarantee sufficient stop.
- Infusion Set: It is convenient to install and uninstall, which ensure precision infusion and safe priming.
- Drainage Bag: closed priming to meet the requirements of quality control, single way drainage bag and double way drainage bag available.
- Customized Designed: different sizes of pumptube and drip chamber to meet the requirements of customers.



Item No.	Size	Cs (PCS)
Hemodialysis Blood Set	Set	1pcs/bag




Respiratory Products

Respiratory Products

Item No.	Size	Cs (PCS)
PAM-001	X-Large , Large ,Medium , S	1pcs/bag
PNM-001	X-Large , Large ,Medium , S	1pcs/bag
POM-001	X-Large , Large ,Medium , S	1pcs/bag
POM-002	X-Large , Large ,Medium , S	1pcs/bag

Constitute



Disposable Anaesthesia Mask

Product Features:

Colored Removable Hook rings for easy size identification and fixation.
Soft inflated cushion for good sealing and patient comfort made of medical grade PVC for patient comfort.



Disposable Nebulizer Mask

Product Features:

- Made of transparent soft PVC for good visualization
- Jet Nebulizer Set come with mask as the interface of inhalation.
- Adjustable nose clip assures comfortable and patient assessment.
- Equipped with a 6cc or 15cc nebulizer Standard Particle size:0.5-1.0 micron
- Size: X-Large (Adult Elongated), Large(Adult Standard) Medium (Pediatric Elongated), S(Pediatric Standard)



Disposable Oxygen Mask

Product Features:

- Made of transparent soft PVC for good visualization
- The turn up rim can ensure comfortable fit with good airtightness.
- Adjustable nose clip assures comfortable and patient assessment.
- Kink-free, crush-resistant tubing
- The standard length of tube is 2.0m,and different length is available.
- Size: X-Large (Adult Elongated), Large(Adult Standard) Medium (Pediatric Elongated), S(Pediatric Standard)



Disposable Oxygen Mask

Product Features:

- Made of transparent soft PVC for good visualization
- The turn up rim can ensure comfortable fit with good airtightness.
- Adjustable nose clip assures comfortable and patient assessment.
- Kink-free, crush-resistant tubing
- The standard length of tube is 2.0m,and different length is available.
- Size: X-Large (Adult Elongated), Large(Adult Standard) Medium (Pediatric Elongated), S(Pediatric Standard)



Endotracheal Tubes

Product Features:

Made from non-toxic PVC, transparent soft and smooth.
Available with cuff(high volume low pressure) X-ray opaque line throughout the length of the tube.
With 15mm standard connector.
Murphy eye at distal end for atraumatic intubation and extubation.
High volume and low pressure cuff provides positive tracheal wall seal.
Designed for use during short and long term ventilation and routine surgical procedures.
Available with Murphy only Packed in individual Blister bag Size: 4.0-10mm



Item No.	Size	Cs (PCS)
PET-001	3mm,4mm,5mm,6mm,7mm,8mm,9mm,10mm	1pcs/bag
PET-002	3mm,4mm,5mm,6mm,7mm,8mm,9mm,10mm	1pcs/bag

Urinary Drainage Bag

Product Features:

Capacity:2000ml
Without bottom outlet; with 90cm inlet tube,
With Anti-Reflux valve,
With pull-push valve,
Sterilization:EO
Packing:1pc/PE bag, 250pcs/ctn
Measurement of carton: 54X 40X27cm



Item No.	Size	Cs (PCS)
Urinary Drainage Bag	54X 40X27cm	1pc/PE bag, 250pcs/ctn

Hydrocolloid Wound Dressing

Product Features:

Hydrocolloid is a material containing gel-form agents, such as carboxymethylcellulose (CMC) and gelatin. It is applied to wound management, support wound bed a moist environment, accelerate wound healing process.



Item No.: 6501
Size: 55×37mm



Item No.: 6502
Size: 30×30mm



Item No.: 6503
Size: 20×65mm



Item No.: 6504
Size: 20×60mm



Item No.: 6505
Size: 42×62mm



Item No.: 6506
Size: 30×30mm



Item No.: 6507
Size: 44×69mm



Item No.: 6508
Size: 80×120mm



Item No.: 6509
Size: 150×180mm



Foam Dressing (Non Adhesive)

Product Features:

The foam dressing Non adhesive is a polyurethane foam with high performance absorption, total exudates management.

Product	Item No.	Size	Package
	6601	5×5cm	20 pcs/box, 10boxes/ctn
	6602	10×10cm	10 pcs/box, 10boxes/ctn
	6603	15×15cm	10 pcs/box, 10boxes/ctn

Foam Dressing (Waterproof ,Non Adhesive)

Product Features:

The foam dressing Waterproof is a polyurethane foam with a highly breathable film backing which prevents wound exudate strike-through and is a barrier to outside contamination while the dressing remains intact without leakage.

Product	Item No.	Size	Package
	6611	5×5cm	20 pcs/box, 10boxes/ctn
	6612	10×10cm	10 pcs/box, 10boxes/ctn
	6613	15×15cm	10 pcs/box, 10boxes/ctn

Foam Dressing (self-adhesive)

Product Features:

The foam dressing Waterproof is a polyurethane foam with a highly breathable film backing which prevents wound exudate strike-through and is a barrier to outside contamination while the dressing remains intact without leakage.

Product	Item No.	Size	Package
	6622	10×10cm	10 pcs/box, 10boxes/ctn
	6633	15×15cm	10 pcs/box, 10boxes/ctn

Film plus

Product Features:

Film Plus is made of polyurethane film, coated with acrylic adhesive with cottonpad, transparent, waterproof, air-permit. Bathing or shower allows while dressing wears.

Product	Item No.	Size	Package
	6201	6×7cm	50 pcs/box, 40boxes/ctn
	6202	10×12cm	50 pcs/box, 24boxes/ctn
	6203	10×15cm	50 pcs/box, 10boxes/ctn
	6204	10×20cm	50 pcs/box, 10boxes/ctn

PORE/Self-Adhesive Wound Dressing

Product Features:

Pore is a soft cloth adhesive wound dressing, hypoallergenic, gentle, used as a primary or second dressing over light to moderately draining wound, securing devices, or for protecting I. V. cannula sites.

Product	Item No.	Size	Package
	6101	7x9cm	50 pcs/box, 20boxes/ctn
	6102	6x8cm	50 pcs/box, 20boxes/ctn
	6103	6x7cm	50 pcs/box, 20boxes/ctn
	6104	10x10cm	25 pcs/box, 20boxes/ctn
	6105	10x15cm	25 pcs/box, 20boxes/ctn
	6106	10x20cm	25 pcs/box, 20boxes/ctn
	6107	10x25cm	25 pcs/box, 10boxes/ctn
	6108	10x30cm	25 pcs/box, 10boxes/ctn
	6109	10x35cm	25 pcs/box, 10boxes/ctn

Strips-Wound skin closures

Product Features:

Strips-wound skin closures are made of non woven, rayon-backed material reinforced with filaments for strength coated with hypoallergenic adhesive. Effectively prevent and reduce the scar proliferate; Sterile package, easy operation and comfort to wear.

Product	Item No.	Size	Package
	6301	3x75mm	5strips/pouch, 50pouches/box, 10boxes/cth
	6302	6x38mm	6strips/pouch, 50pouches/box, 10boxes/cth
	6303	6x75mm	3strips/pouch, 50pouches/box, 10boxes/cth
	6304	6x100mm	10strips/pouch, 50pouches/box, 10boxes/cth

Sports tape - kinesiology tape

Product Features:

Kinesiology tape is made of 98% cotton, 2% nylon fabric coated with acrylic adhesive. Elasticity is more than 1: 1.6, water proof, high adhesive, no glue release, with separate paper. Twelve colors available, can die-cut and print.

Color Bar

White	Blue	Skin	Yellow	Red	Pink	Light Green	Light Blue	Green	Orange	Royal Blue	Black
-------	------	------	--------	-----	------	-------------	------------	-------	--------	------------	-------



Product	Item No.	Size	Package
	6401	25mmx5m	12pcs/box, 20boxes/ctn
	6402	50mmx5m	6pcs/box, 20boxes/ctn
	6403	75mmx5m	4pcs/box, 20boxes/ctn

Silicone Foam Dressing

Product Features:

Silicone Foam is a foam dressing with silicone gel coating, designed to maintain a good healing environment for moderate to high exudate wounds such as pressure ulcers, leg and foot ulcers, traumatic wounds and other secondary healing wounds.

Product	Item No.	Size	Package
	6701	5x5cm	10/box, 400/case
	6702	7.5x7.5cm	10/box, 400/case
	6703	10x10cm	10/box, 400/case
	6704	10x20cm	5/box, 200/case
	6705	10x30cm	5/box, 200/case
	6706	15x15cm	10/box, 400/case
	6707	15x20cm	5/box, 200/case
	6708	17.5x17.5cm	5/box, 200/case
	6709	Heel 12.5x20cm	5/box, 200/case
	6710	M Sacral 18x18cm	5/box, 200/case
	6711	L Sacral 23x23cm	5/box, 200/case



First Aid Bag

First Aid Bag

Product Features:

IFAK (Individual First Aid Kit) is designed for individual soldier self and mutual rescue in the most valuable 10 minutes of battlefield treatment. It has the functions of hemostasis, binding, creating nasopharyngeal airway, auxiliary fix, pain relieving, anti-infection, etc. There is Zeolite Hemostatic Gauze inside the kit, and it is waterproof and flame retardant, easy to carry and use, suitable for variety of locations. It is the necessary First Aid Device in individual combat.



Camouflage



Red



Black



Tan



Green



Outdoors Kit



Car Kit



Office Kit



Family Kit



First Aid Kit

Optional list

Sterilize	Dressing	Binder	Physiotherapy	Medical Supplies	
Bandage	Sterile cotton wool	Tourniquet	Fever Cooling Patch	Tweezers	ARS Needle Kit
Soap	Compressed Gauze	Rolled Splint	Burn Dressing	Nitrile gloves	Pocket mask
Lodine Tablets	Sterility Dressing	S-Rolled Gauze	Burn Gel Dressing	Eye Protection	Thermometer
Saline	Trauma Military Field Dressing	Hemostatic tampon bandage		Trama Scissors	
Swab Rod	Chitosan Hemostatic Granules	Elastic sport bandage		Emergency Blanket	
Alcohol Pad	Zeolite Hemostatic Gauze	self-adhering bandage		Permanent marker	
	Rapid Hemostatic Dressing	Israeli Battle Bandages		Nasopharyngeal Airway	
	Chest Seal Vented	Triangular Bandages Non-Woven		Mouth to Mouth Barrier	
	Kaolin Hemostatic Gauze			Nasopharyngeal airway with lubricat	



Tourniquet



Triangular Bandages Non-Woven



Israeli Battle Bandages



Zeolite Hemostatic Gauze



Kaolin Hemostatic Gauze



Hemostatic Powder



Compressed Gauze



Chest Seal Vented



Mouth to Mouth Barrier



Non-Adherent Pads



Rolled Splint



Emergency Blanket



Hemorrhage Control Bandage



Self-Adhering Bandage



Nitrile gloves




Trauma Scissors



ARS Needle Kit




Permanent marker



Nasopharyngeal Airway

1 Tourniquet



1. Place the wounded limb through the tourniquet band loop.
2. Pull the band very tight, Fasten it back and adhere.
3. Twist the rod until bleeding has stopped, Place the rod inside the clip and lock it in place.
4. Secure the rod and record the time of application on White Strap.



2 Israeli Battle Bandages



1. Open the package and place pad on wound.
2. Wrap the elastic bandage around limb or body part tightly.
3. Insert elastic bandage into pressure bar, Pull back-forcing pressure bar down onto the pad, continue to wrap tightly.
4. Secure hooking ends of closure bar into elastic bandage.



3 Zeolite Hemostatic Gauze



1. Clean the wound.
2. Open package and pack zeolite dressing into wound.
3. Apply pressure directly over bleeding source.
4. Wrap and tie bandage maintain pressure. Seek medical care immediately.



4 Kaolin Hemostatic Gauze



1. Identify the wound and expose the injury.
2. Pack the most active bleeding wound first.
3. Apply direct, firm pressure for at least 3 minutes or until the bleeding stops.
4. Reassess to ensure bleeding has stopped.





Pet First Aid Bag

Product Features:

- Special tapered &kink resistant FEP radiopaque catheter for easy insertion with optimal flow rate
- Special flexible wings to provide easy & secure fixation
- Manufactured from tested Biocompatible materials offering longer indwelling time
- Rounded grip on the needle hub for better support during insertion.
- Injection port with non-return silicon valve for intermittent medication
- Transparent flashback chamber for easy visual medication visual medication



Item No.	Size	Cs (PCS)
Pet First Aid Bag	Large,Medium,Small	1pcs/bag

Pet Medical



Pet oral cleaning nutrient solution

Wound dressing for pets

Item No.	Size	Cs (PCS)	Available Color
Pet oral cleaning nutrient solution	30ml	1pcs/bag	Yellow

Item No.	Size	Cs (PCS)	Available Color
Wound dressing for pets	10ml	1pcs/bag	Blue

Pet First Aid Bag

Pet Medical



PLO-50



PAB-100




PVG-01

Item No.	Size	Cs (PCS)	Available Color
PLO-50	0.5ml	6pcs/box, 24boxes/ctn	Blue
PLO-100	1.0ml	6pcs/box, 24boxes/ctn	Blue

Item No.	Size	Cs (PCS)	Available Color
PAB-050	0.5ml	1pcs/box, 100boxes/ctn	Blue
PAB-100	1.0ml	1pcs/box, 100boxes/ctn	Blue
PAB-200	2.0ml	1pcs/box, 100boxes/ctn	Blue

Item No.	Size	Cs (PCS)	Available Color
PVG-01	3.0ml	1pcs/box, 100boxes/ctn	Blue




Certifications




ISO13485:2016 Certificate



Testing Report



Foreign Trade Record Form



Medical Device Business Record form



Business License



Awarded As A High-tech Enterprise Certificate

Quality Control



Sales Displays



Lined writing area with a vertical blue line on the right side and a series of punch holes on the left side.

NOTES

NOTES:

Any product specifications and recommendations given herein must not be seen as guaranteed product characteristics. They are based on our laboratory tests and on practical experience. Since individual application conditions are beyond our knowledge, control and responsibility, this information is provided without any obligation. We do warrant the continuously high quality of our products being free from defects in accordance with and subject to our General Sales Conditions. However, own adequate laboratory and practical tests to find out if the product in question meets the requested properties are recommended. A claim cannot be derived from them. The user bears the only responsibility for non-appropriate or other than specified applications.